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Endogenic
bio-resonance
correction device

DeVita BRT

Manual • Passport

Purpose of the device

The «Devita BRT» device (from now on referred to as the device) is a generator of weak electromagnetic radiation, operating in the range of natural oscillations of organic systems. It is intended to strengthen its own defensive forces and to correct natural electromagnetic emissions of natural biosystems. The action is based on the method of F.Morrel.

Attention! The device can be used by children only under the supervision of adults. Children are not allowed to play with the device. Individuals with limited physical abilities and mental disabilities may use the device only under the supervision of the appropriate personnel.

Content

Device «DeVita BRT» 1 pc.
Cylindrical manual electrode 2 pc.
Platelet foot electrode 2 pc.
Set of wires 1 pc.
Charger «Li500-2» 1 pc.
Battery Li-Ion (14500) 3.7B\900mAh 2 pc.
The network adapter 220V\50Hz\+7,5V\2A 1 pc.
Instruction manual and passport 1 pc.

Specifications

1. Method F. Morrel
2. Number of leads, pcs. 4
3. Gain, not less than 50
4. Frequency range, kHz 0-180
5. Power supply:

AC adapter 220V / 50Hz \ + 7.5V / 2A
Battery Li-Ion (14500) 3.7B\900mAh

6. Overall dimensions 117x79x33mm
7. Weight, without accessories 0.2kg

Safety Information

Device:
- safe to use
- not toxic
- not flammable
- not explosive
- not a source of radiation hazardous to health.

Environmental Protection

Rechargeable batteries can contain toxic substances that are harmful to health and the environment. Therefore, after the end of their service life, do not throw them into household garbage. Remove the batteries from the device (in accordance with the item «Preparation for work») and transfer it to a specialized point for further disposal. By doing this you can protect the environment.

Discharged batteries can damage the device. If you do not use the unit for a long time, remove the batteries from the battery compartment. When batteries leak out, wear protective gloves and clean the battery compartment with a dry cloth. Dispose of the removed batteries in a specialized place for further disposal.

Rules of care, storage and transportation

The device without packaging should be stored indoors at a temperature of 10 to 35° C and relative humidity of not more than 80%.

The device must be protected from direct exposure to moisture and the action of active chemicals.


Do not disassemble the device.

Do not drop the device.

Do not allow anything to fall on the unit.

We recommend that you keep the packaging during the warranty period.

Device Description

1. Indication of connected leads
2. Status indication («PPROCESS» is the procedure being performed, «LOW BATT» is the battery discharge, «TEST» is the «Test» mode)
3. «TEST» - the button for activating the «Test» mode
4. «ON» - power on / off button
5. «RUN» - the button for activating the «Session» mode
6. Connector for a standard network adapter
7. Connector (LH) for connecting the left cylindrical hand electrode
8. Connector (RH) for connecting the right cylindrical hand electrode

9. Connector (LF) for connecting the left plate foot electrode

10. Connector (RF) for connecting the right plate foot electrode

11. Hand electrodes

12. Foot electrodes

Attention! The battery compartment is located on the back of the device.

Terms of use

1. Do not wipe the surface of the device with alcohol and other organic solvents, remove dirt with light efforts using a soft brush or a cotton pad with water.

2. Do not place the device near extremely warm, dusty or damp places.

3. After storing the device in a cold room or after transporting it to winter conditions, it can be switched on, only after 1 hour of residence at room temperature.

Working with the device

Preparation for work

1. Check the presence of batteries in the battery compartment, for which:

- turn the device upside down
- press your finger on the surface of the battery cover and slide it
- Make sure the batteries are present.

2. Close the battery compartment.

3. When using the network, first connect the device to the network adapter, then connect the adapter to the network. In this mode, the batteries do not work.

4. Connect the hand and foot electrodes to the corresponding connectors of the device using the standard set of wires.

5. Switching on is done by pressing the «ON» button on the device panel.
The green «PROCESS» LED will flash.

6. Turn off the device by pressing the «ON» button. In this case, the «PROCESS» LED should turn off.

7. If necessary, perform a performance check, for which:

- press the «TEST» button. The «TEST» LED and the «LH-RH» lead lights should light up. Connect the LH electrodes together with the RH. There must be a continuous sound;

- by pressing the «TEST» button and connecting the «LH-LF», «LF-RF», «RF-RH»

lead electrodes together, make sure that a continuous sound is heard at each connection and the corresponding LEDs are lit;

- if there is no sound on some lead, check the electrode-cable-device connection;

- If there is no sound after that, then this lead does not work. The device can not be used;

- if necessary, check any lead without having to go through the previous ones, press the «TEST» button and hold it until the LEDs of the desired lead light up, and then release;

- after checking the lead «RF-RH», open all the leads and press the «TEST» button again. There must be a continuous sound. If there is no sound, the device can not be used.

Attention! This procedure should be performed after every 10 sessions.

8. If, during battery operation, the "LOW BATT" LED begins to flash and a sound comes on, you need to turn off the device using the "ON" button, remove the batteries and charge them.

9. The batteries are charged as follows:

- turn the device upside down;

- press your finger on the surface of the battery cover and slide it;

- remove the batteries, install them in the «Li5002» charger and charge according to the description on the «Li500-2»;

- install charged batteries in the battery compartment and close the battery compartment cover.

Conducting a prophylactic session

The device can be used in two modes: «FM (fast mode)» and «BM (base mode)».

Attention! During the session, there should be no direct contact between the arms and the legs. You can have contact through the fabric.

1. FM mode

The session is conducted only on one (any) lead.

1.1. To select a mode, install a short cable to the «RF-LF» lead.

1.2. Connect the electrodes of the selected lead to the «LH-RH» sockets.

1.3. Click the «RUN» button. The "PROCESS" LED should be lit permanently.

Within the next 10 seconds, take the electrodes in your hands and wait for the start-sound and the "LH-RH" LEDs to turn on (the session has started). After 5 minutes, the session will end, the end-sound will declare that the session is over.

If there is no contact in the circuit after 10 seconds after pressing the «RUN» button, for example, the hand-electrode-cable-device (bad contact of the electrode with the hand), the lead LEDs flash, the LED «PROCESS» does not make a 'beep' sound.

If the cause of the open circuit, is resolved in less than 10 seconds, the session will automatically continue. If the cause is not resolved within 10 seconds, the device automatically returns to its original state.

1.4 At the end of the session, turn off the device with the "ON" button or, after 20 seconds, the device will automatically turn off.

2. The BM mode

The session is conducted over all leads.

2.1. This mode is set by default, if there is a short cable in the lead «RF-LF».

2.2. Connect the left and right cylindrical manual electrode to the corresponding sockets «LH» and «RH». Connect the left and right plate-foot electrode to the corresponding «LF» and «RF» sockets and place bare feet on them.

2.3. Press the «RUN» button. The «PROCESS» LED should light up. Within the next 10 seconds, pick up the electrodes and rain the audio signal «trill up» and turn on the lead LEDs «LH-RH» (the session started). The leads «LH-RH», «LH-LF», «LF-RF», «RF-RH» must be switched in this order with an interval of 5 minutes for each lead. After 20 minutes, the session will end, the end-sound will declare that the session is over.

If there is no contact in the circuit 10 seconds after pressing the «RUN» button, for example, the hand-electrode-cable-device (bad contact of the electrode with the handle), the LEDs of the corresponding lead blink, there is no switching to the next lead, the «PROCESS» LED does not light up and a sound signal sounds. If the cause of the chain termination is eliminated in less than 10 seconds, the session will automatically continue. If the cause is not resolved - within 10 seconds the device will automatically reset.

3. At the end of the session, turn off the device with the "ON" button or, after 20 seconds, the device will automatically turn off.

Manufacturer's Warranty

The manufacturer guarantees that the «DeVita BRT» device conforms to the technical specifications when the user observes the rules of operation, transportation and storage.

The warranty period of the device is 24 months from the date of sale through the retail network.

If there is no date of sale and the stamp of the trade organization in the passport, the warranty period is calculated from the date of manufacture of the device by the manufacturer.

During the warranty period, the owner is entitled to free repair. Warranty repair is carried out on the territory of the manufacturer. Delivery of the defective device to the repair site is carried out at the expense of the buyer.

If there is a violation of the safety of the device case and / or in the absence of a warranty card, no claims are accepted and no warranty repair is performed.

The following faults do not apply to the warranty:

- defects due to improper operation;
- defects caused by natural disasters;
- in violation of the integrity of the device sealing;
- in the presence of mechanical damage (cracks, chips, etc.).

The following products are not covered by the warranty:
- charger;
- batteries;
- AC Transformer.

The buyer has the right to replace the defective device with a new one during the warranty period in the following cases:

- when repairing the device for the fourth time;
- the device is not subject to repair.

Attention! The manufacturer reserves the right, without prior notification, to make changes in the design of the device and the software without impairing the existing characteristics.

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EC Declaration of Conformity
Nr.W16

CE

The company **DETA-ELIS EUROPA GMBH** with a seat at Kaiserstr. 74, 63065 Offenbach, Germany

Herewith to DECLARE

That the equipment described below complies with all applicable at the time of construction essential requirements of the **Directive 2011/65/EU (RoHS)** and the **Directive 2014/30/EU (EMC).**
This declaration relates exclusively to the equipment in the state in which it was placed on the market, and excludes alterations or manipulations, components which are added and/or operations carried out subsequently by the final user without our official consent.

Equipment: DeVita
Type / Model: BRT
Standards Applied: IEC 61000-4-2-2008
GOST P МЭК 60601-2:2014
GOST 52770:2007
GOST ISO 10993-1-2-5-10-11-12-18-2011

Date: 01 July 2018

Authorized Representative: The General Manager Ms. Konopleva Tatyana

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Deta-Elis Europa GmbH